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浙江安吉赛安芙生物科技有限公司
2nd Floor, No. 3 Factory, No. 489 Wenyun Road, Tangpu Industrial Park, Dipu Subdistrict, Anji County, Huzhou City, Zhejiang Province, 313303, P.R. China. | www.safheal.com

Reference (clinical background): FDA. FDA Clears First Blood Test Used in Diagnosing Alzheimer's Disease (2025). FDA publication. The FDA reference describes a different assay.

Aβ42 Rapid Test

For professional in vitro use | Human serum / EDTA plasma

Pack sizes: 1, 5, 10, 20 or 25 tests/box; see box label | Specimens: EDTA plasma or serum

Intended Use

For qualitative detection of elevated amyloid-beta 42 (Aβ42) in human EDTA-anticoagulated plasma or serum as an aid to assessment. Interpret the result with the medical history, clinical findings and other appropriate tests. This test alone cannot diagnose or exclude Alzheimer's disease.

Summary

Aβ42 is a biomarker associated with amyloid pathology research. Results depend on the analytical method, specimen type and pre-analytical handling; values from different methods or matrices are not automatically interchangeable. This single-analyte test is qualitative and does not provide an exact Aβ42 concentration.

Test Principle

The test uses a colloidal gold immunochromatographic double-antibody sandwich method. Aβ42 in the specimen binds gold-labeled anti-Aβ42 antibody, migrates along the membrane and is captured by a second anti-Aβ42 antibody coated at T to form a visible band. C is an internal procedural control. With a valid C line, C and T together indicate elevated Aβ42.

Reagents and Materials

Pack sizes: 1, 5, 10, 20 or 25 tests/box; see the box label and packing list for the actual configuration. Each test uses a test device, 1 drop of specimen and 3 drops of matching buffer. Also required: timer, specimen labels and personal protective equipment.

Specimen Collection and Handling

Use EDTA-anticoagulated plasma or serum. Trained personnel should collect venous blood using standard procedures: use an EDTA collection tube for plasma and a serum tube for serum. Separate specimens according to the tube manufacturer's instructions and laboratory procedures; test the separated plasma or serum. Identify each specimen and prevent cross-contamination. Follow product-validated specimen stability, storage and transport requirements. Do not substitute whole blood or plasma collected with another anticoagulant for the confirmed specimen types.

Storage and Stability


Follow the storage conditions and expiry date on the product label. Keep sealed, dry and away from direct sunlight. Do not use expired, damp or damaged products. Open immediately before use and test promptly; do not mix reagents from different lots or products.

Warnings and Precautions

For trained personnel for in vitro testing. Handle human specimens as potentially infectious, use suitable personal protective equipment and prevent cross-contamination. Do not reuse test devices or single-use materials. Use the supplied matching buffer and do not mix reagents from different lots or products. Follow the application procedure and reading time; dispose of waste according to applicable requirements.

Specimen Collection and Test Procedure

01 Collect and separate

Trained personnel collect venous blood and separate per tube and lab instructions. Use only separated EDTA plasma or serum.

02 Add specimen

1 drop

Lay the cassette flat. Add 1 drop of separated EDTA plasma or serum to the single S well.


03 Add buffer

3 drops

Add 3 drops of supplied buffer to the same S well. Avoid bubbles and start timing immediately.

04 Read and record

10 min

Read and record at 10 minutes.
Do not interpret at 15 minutes or later.

Aβ42 Rapid Test

For professional in vitro use | Human serum / EDTA plasma

Pack sizes: 1, 5, 10, 20 or 25 tests/box; see box label | Specimens: EDTA plasma or serum

Test Procedure

1. Check package integrity and expiry. Prepare reagents and specimen according to product instructions, identify the device and specimen, and place the cassette level.
2. Add 1 drop of EDTA-anticoagulated plasma or serum to the only S well, followed by 3 drops of the supplied buffer. Avoid bubbles and start timing immediately after adding the buffer.
3. Observe and record the result at 10 minutes after application. Results at 15 minutes or later are invalid and must not be interpreted. If C is absent, repeat using a new cassette.

Quality Control

C must appear in every valid test; absence of C makes the result invalid. C monitors migration but does not establish adequate specimen collection or correct Aβ42 concentration. Use appropriate external controls according to laboratory quality procedures.

Interpretation of Results

C + T: elevated indicated

C only: elevated not indicated

No C: invalid

C and T present: elevated Aβ42 is indicated under this test's reading rule; read any visible T line, including a faint line. C only, with no T: elevated Aβ42 was not indicated by this test. C absent: invalid regardless of T; repeat using a new cassette. None of these results alone diagnoses or excludes Alzheimer's disease. Do not infer an unvalidated quantitative concentration from band intensity.

Limitations and Testing Recommendations

This test cannot alone diagnose or exclude Alzheimer's disease or independently establish cerebral amyloid pathology. The concentration cutoff for a positive result is not specified; a result cannot establish whether the clinical reference range is exceeded. Interpret with quantitative testing and clinical information. Results may be affected by sampling, pre-analytical handling, methods and patient clinical status; investigate discordant, invalid or difficult-to-read results.

Comparison with a Chemiluminescence Assay

A total of 160 specimens were compared with a chemiluminescence method: 60 comparator-positive and 100 comparator-negative; the comparator's classification threshold remains to be confirmed. The following figures describe agreement for the supplied groups, not diagnostic sensitivity, specificity or accuracy for Alzheimer's disease.

Test result	Comparator positive	Comparator negative	Total
Positive	59	2	61
Negative	1	98	99
Total	60	100	160

Positive percent agreement: 98.33% (95% CI 91.14%-99.71%)

Negative percent agreement: 98.00% (95% CI 93.00%-99.45%)

Overall percent agreement: 98.13% (95% CI 94.63%-99.36%)

Two-sided 95% Wilson confidence intervals; overall percent agreement rounded half-up.

Zhejiang Anji Saianfu Biotech Co., Ltd.
2nd Floor, No. 3 Factory, No. 489 Wenyun Road, Tangpu Industrial Park, Dipu Subdistrict, Anji County, Huzhou City, Zhejiang Province, 313303, P.R. China. | www.safheal.com

Reference (clinical background): FDA. FDA Clears First Blood Test Used in Diagnosing Alzheimer's Disease (2025). FDA publication. The FDA reference describes a different assay.