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Canine/Feline Neuter Status Rapid Test

Gonadal tissue marker | Canine/feline | Serum preferred

Cat.no: 0010141

Colloidal gold immunochromatography | Visual qualitative test

One C/T window · One S well · Whole blood / serum / plasma

Proprietary Technology Statement

To protect the proprietary recognition system and core technical information, this IFU does not disclose the molecular identity of the detection marker or key recognition components. Interpret results within the intended use, reading criteria and limitations stated in this IFU.

Intended Use

For qualitative detection of a gonadal-tissue-associated marker in canine/feline whole blood, serum or plasma, to aid assessment of unknown neuter history, suspected ovarian remnants, cryptorchidism or testicular remnants. Interpret with history, examination and other findings. The test alone does not establish neuter status, complete tissue removal or fertility.

Summary and Clinical Context

Ovarian or testicular tissue may be associated with a detectable circulating marker. Prior surgery, sexual maturity and tissue function affect interpretation. Typical scenarios include males with no scrotal testes and unknown history, and females showing estrous signs after spaying.

Test Principle

The assay uses colloidal gold immunochromatography. The marker reacts with gold-labeled recognition components and migrates to T, where immobilized recognition components capture it to form a visible colored line. C monitors the flow process. The designated positive decision concentration is ≥0.20 ng/mL. Report a discernible T line as detected; do not calculate concentration from its intensity.

Contents and Required Materials

Test cards, matched buffer and IFU; see the packing list for quantities and sampling-device configuration. Also required: a pipette or calibrated transfer device capable of measuring 30 µL accurately, suitable collection/separation supplies, timer and protective equipment.

Storage and Expiry

Store under the product-label conditions; see packaging for lot and expiry. Keep sealed and dry. Open immediately before use. Do not use expired, damp or damaged kits.

Specimen Selection and Handling

Serum is preferred. Apply 30 µL for each specimen type: whole blood, serum or plasma. Record animal identity, sex, age/maturity, neuter history, surgery date and reproductive medication.

Serum: allow clotting, separate by routine procedures and transfer clear supernatant without cells or fibrin. Whole blood: collect into an anticoagulant tube confirmed compatible with this product; mix gently before sampling and reject clotted specimens. Plasma: anticoagulate and separate appropriately, avoiding the cell layer.

Use fresh specimens where possible. Avoid contamination, marked hemolysis or severe lipemia. Investigate clots, obscuring background or abnormal flow and recollect/retest. Delayed testing requires handling, storage and transport confirmed suitable for this product. Do not predilute.

Timing and Testing Recommendations

After surgery: wait at least one month after ovary/testis removal before sampling, to reduce effects of circulating residual marker.

Age: female interpretation primarily concerns sexually mature animals; maturity varies by breed and individual. Interpret young or reproductively senescent animals cautiously; an intact older female may test negative.

Ovarian remnants: predominantly luteal tissue may produce progesterone while the marker remains low. Recurrent estrous signs after spaying warrant progesterone and other appropriate assessment even with a negative test. Consider exogenous progesterone exposure when interpreting an elevated result.

Clinical Interpretation (C Valid)

Result	Suggested report and interpretation limits
Marker positive	Suggests ovarian or testicular tissue is present. Alone, it cannot distinguish intact ovaries from remnants, locate a retained testis or establish fertility.
Marker negative	No signal meeting this test's reading threshold was detected. Under applicable conditions this is consistent with neutering, but does not prove complete removal or exclude remnants.
Difficult to interpret	Check specimen and procedure, then repeat or obtain laboratory quantitative confirmation. Do not force a "neutered" classification. A clearly discernible faint T line remains positive.

Laboratory Quantitative Confirmation Intervals

Concentration (ng/mL)	Confirmation category
<0.10	Negative
0.10 to <0.20	Indeterminate
≥0.20	Positive

These are the designated confirmation intervals, not concentrations measured by this strip. The absence of a T line does not distinguish <0.10 from 0.10 to <0.20. Use a suitable, comparably calibrated method and laboratory interpretation.

Zhejiang Anji Saianfu Biotech Co., Ltd.
2nd Floor, No. 3 Factory, No. 489 Wenyun Road, Tangpu Industrial Park, Dipu Subdistrict, Anji County, Huzhou City, Zhejiang Province, 313303, P.R. China. | www.safheal.com

Canine/Feline Neuter Status Rapid Test

Cat.no: 0010141

Colloidal gold immunochromatography | Visual qualitative test

One C/T window · One S well · Whole blood / serum / plasma

Test Procedure

1. Check specimen identity, packaging and expiry. Bring the card, matched buffer and prepared specimen to room temperature. Prepare the specimen as directed; serum is preferred. Open the pouch and place the card level.
2. Accurately transfer 30 µL of whole blood, serum or plasma to the sole S well. Gently mix whole blood before sampling. Avoid bubbles and keep liquid out of the result window.
3. Then add 2 drops of matched buffer, 80 µL total, to the same S well. Hold the buffer bottle upright without touching the card. Start timing immediately after adding the buffer.
4. Keep the card level; do not move, tilt or add more liquid. At 10 minutes, check C first and then read T in good light. Results after 15 minutes are invalid.

Procedure Overview

Serum preferred; 30 µL for every specimen type. Add specimen, then buffer, to the same S well.


01 Collect and prepare

Serum (preferred) | Whole blood | Plasma

Prepare the appropriate specimen

02 Add specimen

30 µL specimen → S

03 Add buffer

2 drops buffer (80 µL) → same S

04 Time and read

10 min

>15 min


✕

Read at 10 min; invalid after 15 min

Result Illustrations and Reading

C = Control | T = Test

Marker detected / Positive

Both C and T appear. Any clearly discernible T line, including a faint line, is reported as positive.

Marker not detected / Negative

C appears and T is absent. No signal meeting the reading threshold was detected; this does not mean “completely neutered.”

Invalid

No C, whether T appears or not, is invalid. Do not interpret obscuring background or abnormal flow. Investigate and repeat using a new card.

Illustrations are schematic. Review unreadable results; a clearly visible faint T line is positive and is not assigned to the quantitative indeterminate interval because it is faint.

Quality Control

C must appear in every test. It is a procedural control, not proof of adequate sampling or correct specimen concentration. Use appropriate external controls under laboratory policy. Never report invalid results as negative.

Precautions

For trained veterinary personnel; single use. Accurately measure 30 µL; one drop from an arbitrary dropper is not necessarily 30 µL. Use only the matched buffer; do not substitute other reagents. Wear gloves and handle/dispose of specimens and waste as potentially infectious.

Limitations

The test gives no exact concentration. Line intensity cannot estimate tissue amount or location. Timing, near-threshold levels, tissue function and specimen quality may affect results. Discordance with history or examination warrants repeat testing and suitable laboratory/imaging assessment. Do not base a decision for or against surgical exploration solely on this test.

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